



Sample Collection Materials Checklist and Ordering Information

- ❑ **PPE according to your SOP**
- ❑ **Copan swabs (at least 2)**

Fisher Scientific 1-800-766-7000 www.fisherscientific.com Item: Enviro-pbstl swab 916C Catalog # NC9053427	Copan Diagnostics, Inc. Toll free: 800-216-4016 909-549-8793 ext. 208 www.copanusa.com Item: Enviro-pbstl swab 916C
--	---
- ❑ **Plastic scoop or spoon (optional)**
- ❑ **Plastic Vial with Cap (no glass) (optional)**
 - Fisher Scientific
1-800-766-7000
www.fisherscientific.com
 - Item: 50 mL tubes, Catalog number 05-539-6
- ❑ **Gallon-sized sealable bags (Ziploc or similar)**
- ❑ **Gallon-sized paint can (Teflon coating recommended)**
 - Empty cans may be purchased at Home Depot, Lowe's, paint and hardware stores
- ❑ **0.5% Bleach solution/spray bottle**
 - To prepare: mix 1 part household bleach with 9 parts water
 - This must be prepared daily or on-scene to be effective
 - Fisher Scientific
1-800-766-7000
www.fisherscientific.com
 - Item: Dual sprayer with water bottle, bleach cartridge starter kit/refills
 - Catalog numbers: 23-640-127, 23-640-128, 23-640-129
- ❑ **Environmental Sample Submission Form**
 - This can be found at the Department of Health website:
<http://www.doh.state.fl.us/lab/laboratoryservices.htm>
- ❑ **Pen and/or marker**